

Ribbond Wetting Resin!

One of the first steps in using Ribbond is to wet the Ribbond with a plain unfilled enamel bonding resin. Most one bottle bonding resins contain organic solvents, dentin primers or self-etchants; these types of resins should not be used to wet the Ribbond. Use your favorite bonding system to bond to tooth structure but wet the Ribbond with a enamel bonding resin that does not contain solvents, primers and etchants. If your office does not have a suitable enamel bonding resin, try *Ribbond Wetting Resin*.

Ribbond® Wetting Resin

- * No alcohol, no organic solvents, primers, or acidic monomers
- * Unfilled resin for optimized fiber bonding
- * Works with all light and dual-cure resins
- * Can also be used as a instrument lubricant



Only \$28.00

**See our technique videos and more
at www.ribbonbond.com**



Ribbond products are sold directly by Ribbond, Inc. 800-624-4554

051518